

# Concentration procedures

A routine part of complete stool examination for parasites Allows detection of small numbers of parasites that may be missed by using only a direct wet smear

Designed to separate protozoan organisms & helminth eggs & larvae from fecal debris by centrifugation and/or differences in specific gravity

There are two types of concentration procedures:

## Sedimentation & Floatation

### 1-SEDIMENTATION METHOD

#### Background:

- Concentration procedures are used to concentrate eggs, larvae and protozoan cysts to increase the sensitivity for parasite detection in stool samples, specially in cases of mild infection.

All type of eggs and cysts can be recovered by sedimentation. Parasites settle down more rapidly by centrifugation.

Larger food particles can be removed prior to centrifugation by filtering through a sieve with a pore size enough to retain those particles. The efficiency of detection is increased by adding formalin for fixation & preservation of parasites, and ethyl acetate to remove organic material especially fat.

#### Materials:

- 10% formalin
- Ethyl acetate
- Centrifuge tubes
- Centrifuge stand
- Funnel
- Gauze
- Spatula
- Pipettes
- Microscope slides & coverslips



## Method:

Step 1:



Emulsify 0.5-1 g of stool in 7 ml of 10% formalin in a tube.

Step 2:



Pour the stool emulsion onto a double layer of gauze in a funnel and collect in a beaker

Step 4:



Pour the filtrate into a 15 ml centrifuge tube and add 3 ml of ethyl acetate, mix well (~1 min) by hand.

Step 3:



Wash the stool through the gauze using formalin. This washes the parasites through but filters out the larger pieces of debris.


Step 5:



Centrifuge the mixture for 5 minutes at 500 g (~2000rpm).



Four layers are formed in the tube after centrifugation  
(Formol-ethyl acetate concentration method)




**Step 7:**



**Mix the sediment well**


**Step 6:**



**Loosen the fatty plug at the top of the tube with an applicator stick and invert the tube to discard the supernatant.**


**\* Few drops should be kept with the sediment.**

**Step 9:**




**Place one edge of the coverslip on the drop and carefully pull the drop along. When the coverslip is in position gently let it drop.**

**Step 8:**



**Place a drop on the slide.**

**Step 11:**



- Scan the entire coverslip systematically as illustrated using the 10X objective.
- If you suspect cysts or trophozoites, use higher magnification.
- 1/3 of the slide should be scanned using a higher magnification.

**Step 10:**



**If there is a thick and thin area, the cysts and ova will generally be found in the thick portion of the mount.**

## **Concentration Methods**

### **2- Flootation**

#### **Flootation procedure**

Permits the separation of protozoan cysts & oocysts, and certain helminthes eggs & larvae through the use of a liquid with a high specific gravity. Parasites are recovered in the surface film, and the debris remains in the bottom of the tube.

#### **Flootation procedure (Cont.)**

##### **\* Advantage:**

- Cleaner preparation than does the sedimentation procedure.

##### **\* Disadvantages:**

- Some helminth eggs (operculated eggs and/or very dense eggs such as unfertilized Ascaris eggs) do not concentrate well with the flootation method.

- Protozoan cysts become distorted and difficult to identify


#### **Zinc Sulfate Flootation Method**

1. Using a tongue depressor, emulsify 2 to 3 grams of feces (a size of a grape) into 5 ml of ZnSO<sub>4</sub> solution (S.G 1.18) in a dish.
2. Push the emulsion through the strainer into a 15 ml centrifuge tube then fill the tube with zinc sulfate solution.
3. Centrifuge for 2 min at 1500 - 2000 rpm.
4. Using a headed-rod or loop, remove a sample from solution surface & place on a microscope slide.
5. Add a drop of iodine to stain cysts & ova

6-prepare the slid as the perspiration way and don the examination.



# Zinc Sulfate Flootation Method




Step 10:



If there is a thick and thin area, the cysts and ova will generally be found in the thick portion of the mount.

Step 11:



- Scan the entire coverslip systematically as illustrated using the 10X objective.
- If you suspect cysts or trophozoites, use higher magnification.
- 1/3 of the slide should be scanned using a higher magnification.